

Intelligent Breed



AGAPPE

DESIGN ► DEVELOP ► DELIVER

ANOTHER MASTERPIECE OF INNOVATION FROM AGAPPE



MISPA-*i*₂

Specific Protein Analyzer

- ✓ UCS Technology - Sensitivity & Linearity at its best
- ✓ Smart Card Calibration
- ✓ 3 Step User friendly procedure
- ✓ Auto incubation | Auto Rotation | Auto Mixer
- ✓ NGSP Certified HbA1c Reagents
- ✓ 500 patient results storage
- ✓ 30 Min battery backup



MISPA-i₂

Specific Protein Analyzer

Mispa i2 is an Innovation in itself which can perform wide range of test parameter with in stipulated amount of time with better quality and accurate result.

Mispa i2 uses Nephelometric technology for protein determination and supports the clinical management of a variety of disease states, such as cardiovascular risk, rheumatoid diseases, inflammation, kidney disease and several other specific protein assays.



Technical Specification

Principle	: Nephelometry and Photometry	Key board	: Numeric and Functional Keys (21 keys)
Calibration	: Smart Card (Integrated Circuit Card)	Display	: 128 x 64 Graphical Display
Light Source	: Laser LED 650nm	Printer	: Inbuilt Thermal Graphical Printer
Laser Power	: 5 mW	Barcode Interface	: Hand held (optional)
Photo detectors	: Laser detector	PC interface	: RS 232 Port
Temperature Control	: Dry heating	Power Supply	: 85 - 240 VAC 50/60 HZ
Temperature range	: Controlled 37 degrees Celsius	Power Consumption	: 100W
Mixing	: Automatic by motor		
Zero setting	: Automatic		

Parameters

Product Name	Assay Range
ASO	50 - 1000 IU/mL
CRP	0.5 - 320 mg/L
RF	10 - 100 IU/mL
HbA1c	3 - 13 %
MICROALBUMIN	5 - 200 mg/L
Lp(a)	1 - 100 mg/dL
CYSTATIN C	0 - 10 mg/L
FERRITIN	1 - 1000 ng/mL
CERULOPLASMIN	5 - 100 mg/dL
ApoA1	15 - 300 mg/dL
ApoB	20 - 330 mg/dL
C3	30 - 400 mg/dL
C4	3 - 80 mg/dL

Product Name	Assay Range
D-Dimer	0 - 10 ug/mL
IgG	360 - 3000 mg/dL
IgA	70 - 630 mg/dL
IgE	5 - 1000 IU/mL
IgM	15 - 400 mg/dL
PREALBUMIN	9 - 75 mg/dL
ALPHA 1 ACID GLYCOPROTEIN	10 - 250 mg/dL
ALPHA 1 ANTI TRYPSIN	40 - 400 mg/dL
ALPHA 2 MACROGLOBULIN	95 - 760 mg/dL
HAPTOGLOBIN	15 - 350 mg/dL
KAPPA LIGHT CHAIN	89 - 800 mg/mL
LAMBDA LIGHT CHAIN	40 - 400 mg/dL

Free Accessories

With MISPA-i2

R1 Pipette 50 - 200 µL	1
Sensor Pipette 20 - 200 µL	1
Micro Pipette 2 - 20 µL	1
Cuvette Stand	1
Pipette Stand	1
Printer Rolls	2

With Reagent Kit

SMART CARD	1
Reaction Cuvettes	15 / 30
Micro Tips	45 / 90



Smart Card
Calibration



Certified HbA1c
Reagents



PT. MULTIMEDILAB KARYAMANDIRI
MEDICAL & LABORATORY EQUIPMENTS

GP Plaza Lantai 3 Unit 6 Jl. Gelora II No. 01 RT. 01/01,
Kelurahan Gelora Kecamatan Tanah Abang, Jakarta Pusat
10270 - INDONESIA
Telp. (021) 3972-9191 (hunting). Fax. (021) 3972-9195

AGAPPE

AGAPPE DIAGNOSTICS SWITZERLAND GmbH
Knonauerstrasse 54 - 6330 Cham - Switzerland Tel: +41 41 780 60 10 Fax: +41 41 780 60 11
info@agappeswiss.com | www.agappeswiss.com