

UA-600/600T

Urine Analyzer

Technical Specifications

Throughput	600 tests per hour
Strips	11 items
Test items	Urobilinogen, Bilirubin, Ketone, Blood, Protein, Nitrite, Leukocytes, Glucose, Specific Gravity, Ascorbic Acid, pH
Wavelength	Red (R), Green (G) and Blue (B)
Principle	Reflectance photometer
Display	Large LCD display, Resolution: 240×64
Input	Touch screen (5.7", UA-600T); Keypad (5.3", UA-600)
Memory	2000 results
Printer	Built-in thermal printer, paper width: 57mm, interface for external printer
Interface	Standard RS-232C
Working environment	Temperature 10°C~30°C Humidity ≤85%
Power supply	AC 100~240V~, 50/60Hz
Power	45VA
Dimension	390mm(H)×330mm(W)×210mm(D)
Weight	5.5kg



U-11 Urine Strips

Urobilinogen, Bilirubin, Ketone, Blood, Protein, Nitrite, Leukocytes, Glucose, Specific Gravity, Ascorbic Acid, pH

- Accurate Results
- Immediate Diagnosis
- Convenient
- High sensitivity and specificity
- Result in 6 seconds
- Easy operation



PT. MULTIMEDILAB KARYAMANDIRI
MEDICAL & LABORATORY EQUIPMENTS
GP Plaza Lantai 3 Unit 6 Jl. Gelora II No. 01 RT. 01/01,
Kelurahan Gelora Kecamatan Tanah Abang, Jakarta Pusat
10270 – INDONESIA
Telp. (021) 3972-9191 (hunting). Fax. (021) 3972-9195



Mindray is listed on the NYSE under the symbol "MR"
Mindray Building, Keji 12th Road South, High-tech Industrial Park,
Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 81888998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com Website: www.mindray.com

UA-600/600T

Urine Analyzer



- 600 tests per hour
- Cold light source with high luminosity, long shelf life
- Innovative light sensor, high sensitivity
- Automatic detection of strip placing and testing
- Large LCD touch screen or keypad configuration
- Automatic used strip collection
- Standard RS-232 interface for connection with external computer
- Built-in thermal printer and interface for external printer
- Storage of 2000 test results
- Multifunctional software



mindray
healthcare within reach