

Time-Resolved Immunofluorescence Analyzer Model: D10 (Portable Type)


Quantitative

Reliable

Accurate

Easy

Fast: Incubating 3-15 min




D10 Time-Resolved Immunofluorescence Quantitative Analyzer

DESCRIPTION

D10 is a small, innovative and accurate POCT Immunofluorescence Quantitative Analyzer that has obtained CE IVDR certificate. Combined with magnetic card auto-calibration, it can test different sample types including whole blood, serum and plasma and quickly provide reliable quantitative results.

Suitable for medical and health institution, clinics and doctor's offices, hospital departments: ICU & CCU, laboratory and scientific research, etc.

Full coverage 60+ Parameters: Exclusive Therapeutic Drug Monitoring (TDM), Tumor, Thyroid, Sex Hormones, Diabetes, Bone Metabolism, Cardiac, Heavy Metal, Inflammation, Immunity, Cerebrovascular, Nephritic, Pneumonic, Gastric, Infectious Disease and Others.



CHARACTERISTICS—Quantitative, 60+ Parameters, Single Channel

Accuracy	$R \geq 0.990$
Specimen	serum, plasma, whole blood
Method	Fluorescence Immunochemical Assay
Screen	7-inch color touch screen
Dimension(L*W*H)	265*230*140mm
Weight	2.0 KGS
Printing	Built-in embedded thermal printer
Test Mode	Quick Test Mode and Standard Test Mode
Result Time	To get results in 10s under Quick Test Mode
Memory	100000 tests
Communication Interface	Serial, Ethernet, USB, Mini USB, Micro SD card

Fluorescence Immunoassay POCT Determination Kits—Available in D10

Therapeutic Drug Monitoring (TDM) Test Items		Category	Test Items
Immunosuppressant	Tacrolimus	Rheumatology	Rheumatoid Factor
	Mycophenolic Acid		Anti-Streptococcal Hemolysin O Test anti cycliccitrullinated peptide antibody
Antiepileptic Drugs (AED)	Carbamazepine	Gastric	Pepsinogen I
	Valproic Acid		Pepsinogen II
	Phenytoin	Sex Hormones	Follicle Stimulating Hormone
	Total Risperidone		Luteinizing Hormone
Antibiotic	Vancomycin		Prolactin
	Teicoplanin		Progesterone
	Meropenem		β-HCG
	Linezolid		Anti-Mullerian Hormone
Antineoplasticdrug	Methotrexate		Estradiol
	Paclitaxel		Testosterone
	5-Fluorouracil	Heavy Metal	Cadmium
Cardiac Glycoside	Digoxin		Plumbum
Antiasthmatic Drugs	Theophylline	Immunity	Rheumatoid Factor
Antituberculous Drugs	Rifampicin		Immunoglobulin E
	Isoniazid	Cerebrovascular	S100-β
Antifungal	Voriconazole	Bone Metabolism	25-hydroxy Vitamin D
Biopharmaceuticals	Adalimumab	Vitamin	Vitamin B12
	Anti-Adalimumab	Pneumonic	Chlamydia Pneumoniae IgG Antibody
	Infliximab		Chlamydia Pneumoniae IgM Antibody
	Anti-Infliximab		Chlamydia Pneumoniae/ Mycoplasma Pneumoniae IgM Antibody
Tumor	Total Prostate Specific Antigen	Cardiac	Creatine Kinase Isoenzymes MB
	Free Prostate Specific Antigen		Troponin I
	Alpha-fetoprotein		D-Dimer
	Cancer Antigen 125		H-Fatty Acid Binding Protein
	Carcinoembryonic Antigen		Myoglobin
	Carbohydrate Antigen 153		N-terminal pro-B-type Natriuretic Peptide
	Cancer Antigen 199		Lipoprotein Phospholipase A2
	Cytokeratin 19 Fragment 21-1		Troponin I/Creatine Kinase IsoenzymesMB/Myoglobin
Thyroid	Total Triiodothyronine	Nephritic	Cystatin C
	Total Thyroxine		Microalbumin
	Thyrotropin		Neutrophil Gelatinase-associated Lipocalin
	Free Triiodothyronine		
	Free Thyroxine		


Items in Red color means no CE

Items in Red color means no CE


Category	Test Items	Category	Test Items
Infectious Disease	2019-nCoV IgG/IgM Antibody	Infectious Disease	Dengue NS1 Antigen
	Influenza A/B Viral Antigen		Dengue IgG/IgM Antigen
	2019-nCoV Antigen		Chlamydia Trachomatis IgG
	2019-nCoV Neutralizing Antibody		Chlamydia Trachomatis IgM
	Human Immunodeficiency Virus (HIV 1+2)		Human T-lymphotropic virus
	Hepatitis C Virus(HCV)	Inflammation	C Reactive Protein
	2019-nCoV/ Flu A+B/RSV combo		Procalcitonin
	Helicobacter Pylori (HP)		Serum Amyloid A
	HBsAg	Diabetes	Interleukin 6
	Flu A+B	Others	Glycated Hemoglobin
	Syphilis		Cholyglycine
	salmonella		Calprotectin
			FOB Determination Kit
			INSULIN
			Ferritin

Diagreat Strong R&D Team specializes in developing Exclusive products and keeps developing new products every year. Diagreat always attach importance to research and development capabilities to meet customers' needs for new products.


How To Use---Easy Operation Guide



Step 1 Read the manual




Step 2 Preparation of D10 & Incubator




Step 3 Reading card



Step 4 Sampling



Step 5 Loading



Step 6 Incubation



Step 7 Testing

E-mail: info@diagreat.com
export@diagreat.com

Whatsapp/Mob:
+86-01069502385
+86-18500563281

[Http://www.diagreat.com](http://www.diagreat.com)



P.T. MULTIMEDILAB KARYAMANDIRI
MEDICAL & LABORATORY EQUIPMENTS

GP PLAZA Lantai 3 Unit 6, Jalan Gelora II No. 01 Rt.01/01
Kelurahan Gelora, Kecamatan Tanah Abang, Jakarta Pusat 10270 – Indonesia
Telp. (021)39729191 (Hunting). Fax : (021)39729195 Email : info@multimedilab.com
<http://www.multimedilab.com>

Add: Floor 1-3, Building 19, No. 2 Yard, Huanke
Middle Road, Tongzhou District, Beijing, China.