




- Bench top, Floor Model optional
- Throughput: Random access, 120T/H (constant speed)
- Probe, sample & reagent position: 2 probes, 40 sample positions, 40 reagent positions and 60 cuvettes
- Cooling system: On board cooling by ON/OFF switch
- Washing and optical system: 8 steps washing and 9 wavelengths





- Bench top, Floor Model optional
- Throughput: Random access, 120T/H (constant speed)
- Probe, sample & reagent position: 2 probes, 40 sample positions, 40 reagent positions and 60 cuvettes
- Cooling system: On board cooling by ON/OFF switch
- Washing and optical system: 8 steps washing and 9 wavelengths



PT. MULTIMEDILAB KARYAMANDIRI
MEDICAL & LABORATORY EQUIPMENTS

GP Plaza Lantai 3 Unit 6 Jl. Gelora II No. 01 Rt. 01/01
Kelurahan Gelora Kecamatan Tanah Abang
Jakarta Pusat 10270 - INDONESIA

Telp. (021) 3972-9191 (hunting) Fax. (021) 3972-9195
Website : www.multimedilab.com Email : info@multimedilab.com

